English

Indication

NanoBone® QD is for the reconstruction of bone defects which are not essential for the stability of the bone structure and for augmentation. NanoBone® QD can be used for all fields of application for which autologous cancellous bone could be used.

Properties

NanoBone® QD consists of NanoBone® | granulate S39 (0.6-2.0 mm) and a water soluble carrier material **containing poloxamer**. In NanoBone® | granulate S39, nanocrystalline hydroxyl apatite (HA) is embedded in a silica gel matrix (amorphous SiO₂). The porosity of the granulate is approx. 80 %.

The transmission of infections or diseases of animal or human origin is ruled out as a result of it being manufactured entirely synthetically. NanoBone® QD does not harden after implantation.

Application

NanoBone® QD may only be used by qualified professionals with experience in the required surgical techniques.

The principles of sterile work and patient medication valid for surgical procedures must be observed.

Preparation:

Fragments of bone and soft tissue in the defect must be removed from the bone defect to be treated.

Depending on the nature and location of the bone defect, additional osteosynthetic measures to stabilise the bone grafting material may be required. Assemble the applicator as shown in Figure 1 “Preparing the applicator”.

Implantation:

NanoBone® QD can be used exclusively or in any mixing ratio with autologous cancellous bone.

The area where NanoBone® QD is implanted must be properly secured with rigid fixation. Maximize the contact between the patient’s bone and NanoBone® QD to ensure bone regeneration. NanoBone® does not set in-situ following implantation.

The bone defect is to be completely filled with NanoBone® QD. Avoid over filling as well as compaction of the bone grafting material in the bone defect.

The filled bone defect must be closed with a primary wound closure ensuring that the defect is covered by the periosteum.

Information

For successful bone regeneration and in order to avoid the dislocation of NanoBone® QD, a suitable bone position (cavity) must be provided, and the immobilisation of the bone grafting material is necessary. Micro movements can impair bone regeneration and must be avoided using suitable surgical techniques (e.g. using membranes).

As a result of its high porosity, NanoBone® bone grafting material can look like spongy bone in x-rays and therefore immediately after implantation show a corresponding x-ray contrast.

NanoBone® QD may no longer be used after its use by date has expired.

Warnings

NanoBone® QD is not suitable for load-bearing defects.

There is no data available on the use of NanoBone® QD with pregnant women and children.

To avoid contamination, the sterile bone grafting material may only be removed from the packaging immediately before use and is to be applied in a sterile working environment. Preparation of unused NanoBone® QD is not possible as the destruction of the nanostructure or changes to properties cannot be ruled out during preparation.

Single use only. Do not attempt to reuse and/ or re-sterilize unused NanoBone® QD.

Interactions

Interactions between NanoBone® QD and other medicines and drugs are not known.

Side-effects

Side-effects as a result of using NanoBone® QD are not known.

Contraindications

- Acute and chronic infections in the area of the surgery (soft tissue infections, inflammatory, bacterial bone diseases, osteomyelitis), with antibiotic treatment the user is to decide on the use of NanoBone® QD by means of a risk-benefits assessment

- Serious metabolic diseases such as e.g. serious, uncontrollable or poorly controlled diabetes mellitus

- Calcium metabolism disorders

- Treatment with steroids and other pharmaceuticals which interfere with calcium metabolism

- Immunosuppressive treatment

- Malignant tumors

- Endocrinological bone diseases

- Application in areas of open epiphyseal plates

- Application not in mechanically instable defects (after stabilization of the defect, e.g. with an osteosynthesis plate the application of NanoBone® is possible),

- In the radiated area and/or before a planned radiation

- Bone defects with no or poor vascularisation

Type and content of the packaging

NanoBone® QD is provided in a sterile applicator to meet the demands of the application. This is packaged in a sterile aluminium protection bag (peel-off packaging). This inner bag is packed in an aluminium protection bag (peel-off packaging). NanoBone® QD must no longer be stored once the inner aluminium protection bag has been opened.

Sterility

NanoBone® QD is sterilised using gamma radiation. Before using NanoBone® QD, its expiration date and the impeccable condition of the sterile packaging must be checked.

Glossary

 This product may not be used if the sterile packaging has been damaged

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