NanoBone® | granulate NanoBone® | block

Synthetic, biodegradable bone grafting material

Indication

NanoBone® bone grafting material is for the reconstruction of bone defects and for augmentation.

Properties

NanoBone® bone grafting material is made up of nanocrystalline hydroxyapatite (HA) which is embedded in a silica gel matrix (amorphous SiO₂). To control the speed of biodegradation, various products with different HA/SiO₂ ratios are offered (retardation of biodegradation with reduction in proportion of SiO₂). The transmission of infections or diseases of animal origin is ruled out as a result of it being manufactured entirely synthetically. NanoBone® products are available as granulate or block, to meet the demands of the application. The porosity of the bone grafting material is approx. 80 %.

Application

NanoBone® | granulate and NanoBone® | block may only be used by qualified professionals with experience in the required surgical techniques.

The principles of sterile work and patient medication valid for surgical procedures must be observed.

Preparation:

Fragments of bone and soft tissue in the defect must be removed from the bone defect to be treated.

Depending on the nature and location of the bone defect, additional osteosynthetic measures to stabilise the bone grafting material are required.

NanoBone® bone grafting material may not be implanted into the bone defect dry. For this it is recommended that the granulate or block is mixed or saturated with the patient's own venous blood before implantation. Mixing or saturating with physiological saline solution is also possible.

When contouring the NanoBone® | block and/or bone bedding, the greatest possible contact surface between the bone grafting material and bone must be produced.

Implantation:

The area where NanoBone® is implanted must be properly secured with rigid fixation.

Maximize the contact between the patient's bone and NanoBone® to ensure bone regeneration.

The bone defect is to be completely filled with NanoBone® | granulate. Avoid overfilling as well as compaction of the bone grafting material in the bone defect.

When using the NanoBone® | block, contact between the bone bedding and bone grafting material must be guaranteed using additional measures. One option is to secure the block using bone screws and bone plates. The block is pre-drilled for this purpose. The filled bone defect must be closed with a primary wound closure ensuring that the defect is covered by the periosteum.

Information

Immobilization of the bone grafting material is necessary for successful bone regeneration. Micro movements can impair bone regeneration and must be avoided using suitable surgical techniques.

As a result of its high porosity, NanoBone® bone grafting material can look like spongy bone in x-rays and therefore immediately after implantation show a corresponding x-ray contrast.

NanoBone® bone grafting material may no longer be used after its use-by date has expired.

Warnings

There is no data available on the use of NanoBone® bone grafting material with pregnant women and children.

To avoid contamination, the sterile bone grafting material may only be removed from the packaging immediately before use and is to be applied in a sterile working environment.

Preparation and resterilisation of unused NanoBone® bone grafting material is not possible, as the destruction of the nanostructure or changes to properties cannot be ruled out during preparation.

Single use only. Do not attempt to reuse and/or re-sterilize unused NanoBone®.

Interactions

Interactions between NanoBone® bone grafting material and other medicines and drugs are not known.

Side-effects

Side-effects as a result of using NanoBone® bone grafting material are not known.

Contraindications

Acute and chronic infections in the area of the surgery (soft tissue infections; inflammatory, bacterial bone disease; osteomyelitis); with antibiotic treatment the user is to decide on the use of NanoBone® | granulate by means of a risk-benefits assessment

- Uncontrolled metabolic diseases (e.g. diabetes mellitus)
- Systemic diseases that negatively impact bone regeneration
- Use of medication known to negatively impact bone regeneration
- Immunosuppressive treatment
- Malignant tumours
- Application in areas of open epiphysial plates
- Application not in mechanically instable defects (after stabilization of the defect, e.g. with an osteosynthesis plate the application of NanoBone® is possible),
- In the radiated area and/or before a planned radiation
- Bone defects with no or poor vascularization

Sterility

The bone grafting material is sterilised by irradiation.

The expiry date and integrity of the sterile packaging must be checked before preparation.

A copy of the summary of safety and clinical performance (SSCP) can be requested at infogermany@biocomposites.com.

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Emergo Australia

Level 20, Tower II, Darling Park

201 Sussex Street, Sydney NSW 2000 Australia

NanoBone® | granulate NanoBone® | block

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Biocomposites GmbH