
NanoBone[®] SBX Putty

Instructions for Use

100109.4

Information as of 2025-06-18

CONTENTS

(US) UNITED STATES ONLY	3
(EN) ENGLISH	5
(FR) FRANÇAIS	6
(DE) DEUTSCH	7
(GR) ΕΛΛΗΝΙΚΑ	8
(IT) ITALIANO	9
(PT) PORTUGUÊS	10
(ES) ESPAÑOL	11
(SE) SVENSKA	12
(CZ) ČESKY	13
(DK) DANSK	14
(SK) SLOVENSKY	15
(BG) БЪЛГАРСКИ	16
(EE) EESTI KEEL	17
(HU) MAGYAR	18
(LT) LIETUVIŠKAI	19
(RO) ROMÂNĂ	20
(SI) SLOVENŠČINA	21
(TR) TÜRKÇE	22
(NO) NORSK	23

NanoBone® SBX Putty

Bone Graft Substitute

Product Description

NanoBone® SBX Putty contains NanoBone®, a synthetic, biodegradable bone grafting material that is composed of non-sintered nanocrystalline hydroxyapatite (HA) which is embedded in a silica gel matrix (amorphous SiO₂). The ratio of HA / SiO₂ has been chosen to optimize the rate of biodegradation for treatment of osseous defects in human bone.

NanoBone® SBX Putty does not contain any components of animal or human origin, eliminating the possibility of transmission of infection or disease. NanoBone® SBX Putty is supplied in a sterile polycarbonate applicator and contains NanoBone® granules (0.6-2.0 mm) mixed with a silica gel matrix. NanoBone® SBX Putty does not set in-situ following implantation.

Indications for Use

NanoBone® SBX Putty is an implant intended to fill bony voids or gaps of the skeletal system (i.e., extremities, posterolateral spine, intervertebral disc space, and pelvis). These osseous defects are surgically created or the result of traumatic injury to the bone and are not intrinsic to the stability of the bony structure. NanoBone® SBX Putty resorbs and is replaced with bone during the healing process. When used in intervertebral body fusion procedures, NanoBone® SBX Putty must be used with an intervertebral body fusion device cleared by FDA for use with a bone void filler.

Contraindications

NanoBone® SBX Putty is not designed or sold for any use except as indicated. Do not use NanoBone® SBX Putty in the presence of any contraindication. NanoBone® SBX Putty is contraindicated where the device is intended as structural support in the skeletal system. Other conditions for which NanoBone® SBX Putty is contraindicated include:

- Severe vascular or neurological disease
- Uncontrolled diabetes
- Severe degenerative disease
- Severely impaired renal function
- Hypercalcemia, abnormal calcium metabolism
- Existing acute or chronic infections, especially at the site of the operation
- Inflammatory bone disease such as osteomyelitis
- Malignant tumors
- Uncooperative patients who cannot or will not follow post-operative instruction, including individuals who abuse drugs and/or alcohol

Application

NanoBone® SBX Putty may only be used by qualified professionals with experience in the required surgical techniques. The principles of sterile work and patient medication valid for surgical procedures must be observed. NanoBone® SBX Putty is not intended for load-bearing applications. The area where NanoBone® is implanted must be properly secured with rigid fixation. Maximize the contact between the patient's bone and NanoBone® SBX Putty to ensure bone regeneration. NanoBone® SBX Putty does not set in-situ following implantation. The bone defect must be completely filled with NanoBone® SBX Putty. Avoid overfilling as well as compaction of the bone grafting material in the bone defect. The filled bone defect should be closed with a primary wound closure ensuring that the defect is covered by the periosteum. NanoBone® SBX Putty may appear radiopaque in x-rays, similar to cancellous bone. Do not overfill or attempt to pressurize the bony defect site, as this may lead to extrusion of the product beyond the site of its intended application and damage to the surrounding tissues, or may lead to fat embolization or embolization of the product into the bloodstream.

Instructions for Use

1. Open outer (non-sterile) and inner (sterile) packaging.
2. NanoBone® SBX Putty is designed to be used alone. When used in intervertebral body fusion procedures, NanoBone® SBX Putty must be used with an intervertebral body fusion device cleared by FDA for use with a bone void filler.
3. Implant NanoBone® SBX Putty. Secure the surgical site after implantation to prevent micromotion and/or implant migration. If NanoBone® SBX Putty material is not positioned properly, remove the material and begin again with a new package of material.

Warnings

NanoBone® SBX Putty is not intended for load-bearing uses. There is no data available on the use of NanoBone® SBX Putty with pregnant women or children. Single use only. Do not attempt to reuse and/or resterilize unused NanoBone® SBX Putty.

Adverse Events

The following potential adverse events (AEs) have been reported with the use of FDA-cleared resorbable bone void filler devices:

- Revisions and/or removals
- Superficial wound or deep wound infection
- Pain/discomfort, swelling, redness, fever, inflammation
- Fluid accumulation, wound dehiscence, drainage
- Debridement/irrigation

- Delayed or nonunion, lack of osseointegration, impaired healing, inadequate bone formation
- Altered handling characteristics leading to failure
- Protrusion, dislodgement, migration, or extravasation (leakage)
- Decreased range of motion, loss of motor function, sensory deficit
- Allergic/immune response
- Blood pressure change
- Hematoma
- Cyst
- Death

Sterility

Contents are STERILE by exposure to gamma radiation unless opened or damaged. Do not use if sterile packaging has been damaged. Do not use if expiration date has been exceeded. Do not resterilize.

Storage, Shelf Life and Disposal

Product should be stored between 5-32°C (40-90°F). The expiration date is printed on the label. DO NOT USE NanoBone® SBX Putty AFTER THE EXPIRATION DATE. NanoBone® SBX Putty is environmentally friendly. No special disposal is necessary.

Caution

U.S. Federal Law restricts this device to sale by or on the order of a physician or hospital. Responsibility for proper selection of patients, for adequate training, for experience in the selection and placement of NanoBone® SBX Putty and for the selection of post-operative protocols rests entirely with the physician.

Glossary

Symbol	Description	Source
	Manufactured by.	AAMI/ANSI/ISO 15223-1: 2016, section 5, reference no. 5.1.1
	Contents are STERILE by exposure to gamma radiation.	AAMI/ANSI/ISO 15223-1: 2016, section 5, reference no. 5.2.4
	Do not use if package is damaged.	AAMI/ANSI/ISO 15223-1: 2016, section 5, reference no. 5.2.8
	Do not resterilize.	AAMI/ANSI/ISO 15223-1: 2016, section 5, reference no. 5.2.6
	Single use. Do not reuse.	AAMI/ANSI/ISO 15223-1: 2016, section 5, reference no. 5.4.2
	Store between 40-90 °F.	AAMI/ANSI/ISO 15223-1: 2016, section 5, reference no. 5.3.7
	Consult instructions for use.	AAMI/ANSI/ISO 15223-1: 2016, section 5, reference no. 5.4.3
	CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician.	21 U.S. Code § 353, paragraph (b)(4)(A)

NanoBone® SBX Putty

CE 0482**Synthetic, biodegradable bone grafting material****Indication**

NanoBone® SBX Putty is for the reconstruction of bone defects which are not essential for the stability of the bone structure and for augmentation. NanoBone® SBX Putty can be used for all fields of application for which autologous cancellous bone could be used.

Properties

NanoBone® SBX Putty consists of NanoBone® | granulate S39 (0.6-2.0 mm) and a water soluble carrier material containing poloxamer. In NanoBone® | granulate S39, nanocrystalline hydroxyl apatite (HA) is embedded in a silica gel matrix (amorphous SiO₂). The porosity of the granulate is approx. 80 %. The transmission of infections or diseases of animal or human origin is ruled out as a result of it being manufactured entirely synthetically. NanoBone® SBX Putty does not harden after implantation.

Application

NanoBone® SBX Putty may only be used by qualified professionals with experience in the required surgical techniques. The principles of sterile work and patient medication valid for surgical procedures must be observed.

Preparation:

Fragments of bone and soft tissue in the defect must be removed from the bone defect to be treated. Depending on the nature and location of the bone defect, additional osteosynthetic measures to stabilise the bone grafting material may be required.

Implantation:

- NanoBone® SBX Putty can be used exclusively or in any mixing ratio with autologous cancellous bone.
- The area where NanoBone® SBX Putty is implanted must be properly secured with rigid fixation.
- Maximize the contact between the patient's bone and NanoBone® SBX Putty to ensure bone regeneration.
- NanoBone® does not set in-situ following implantation.
- The bone defect is to be completely filled with NanoBone® SBX Putty. Avoid overfilling as well as compaction of the bone grafting material in the bone defect.
- The filled bone defect must be closed with a primary wound closure ensuring that the defect is covered by the periosteum.

Information

For successful bone regeneration and in order to avoid the dislocation of NanoBone® SBX Putty, a suitable bone position (cavity) must be provided, and the immobilisation of the bone grafting material is necessary. Micro movements can impair bone regeneration and must be avoided using suitable surgical techniques (e.g. using membranes). As a result of its high porosity, NanoBone® bone grafting material can look like spongy bone in x-rays and therefore immediately after implantation show a corresponding x-ray contrast. NanoBone® SBX Putty may no longer be used after its use-by date has expired.

Warnings

- NanoBone® SBX Putty is not suitable for load-bearing defects.
- There is no data available on the use of NanoBone® SBX Putty with pregnant women and children.
- To avoid contamination, the sterile bone grafting material may only be removed from the packaging immediately before use and is to be applied in a sterile working environment.
- Preparation of unused NanoBone® SBX Putty is not possible as the destruction of the nanostructure or changes to properties cannot be ruled out during preparation.
- Single use only. Do not attempt to reuse and/or re-sterilize unused NanoBone® SBX Putty.

Interactions

Interactions between NanoBone® SBX Putty and other medicines and drugs are not known.

Side-effects

Side-effects as a result of using NanoBone® SBX Putty are not known.

Contraindications

- Acute and chronic infections in the area of the surgery (soft tissue infections; inflammatory, bacterial bone disease; osteomyelitis); with antibiotic treatment the user is to decide on the use of NanoBone® SBX Putty by means of a risk-benefits assessment
- Uncontrolled metabolic diseases (e.g. diabetes mellitus)
- Systemic diseases that negatively impact bone regeneration
- Use of medication known to negatively impact bone regeneration
- Immunosuppressive treatment
- Malignant tumours
- Application in areas of open epiphyseal plates
- Application not in mechanically instable defects (after stabilization of the defect, e.g. with an osteosynthesis plate, the application of NanoBone® is possible),
- In the radiated area and/or before a planned radiation
- Bone defects with no or poor vascularisation

Type and content of the packaging

NanoBone® SBX Putty is provided in a sterile applicator to meet the demands of the application. This is packaged in a sterile aluminium protection bag (peel-off packaging). This inner bag is packed in an aluminium protection bag (peel-off packaging). NanoBone® SBX Putty must no longer be stored once the inner aluminium protection bag has been opened.

Sterility

NanoBone® SBX Putty is sterilised using gamma radiation. Before using NanoBone® SBX Putty, its expiration date and the impeccable condition of the sterile packaging and the aluminium protection bag must be checked.

A copy of the summary of safety and clinical performance (SSCP) can be requested at info@germany@biocomposites.com.

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Biocomposites GmbH
Biocomposites Ltd
Biocomposites Inc
My Biologics PTY LTD