NanoBone[®] I granulate

Instructions for Use

100012.8

Information as of 2025-12-04

CONTENTS

(EN) ENGLISH	3
(AU) AUSTRALIA - ORTHOPAEDIC	4
(AU) AUSTRALIA - DENTAL	5
(FR) FRANÇAIS	6
(DE) DEUTSCH	7
(GR) ΕΛΛΗΝΙΚΑ	8
(IT) ITALIANO	9
(PT) PORTUGUÊS	10
(ES) ESPAÑOL	11
(SE) SVENSKA	12
(CZ) ČESKY	13
(DK) DANSK	14
(SK) SLOVENSKY	15
(BG) БЪЛГАРСКИ	16
(EE) EESTI KEEL	17
(HU) MAGYAR	18
(LT) LIETUVIŠKAI	19
(RO) ROMÂNĂ	20
(SI) SLOVENŠČINA	21
(TR) TÜRKÇE	22
(NO) NORSK	23

NanoBone® I granulate

C € 0482

NanoBone® | block
Synthetic, biodegradable bone grafting material

Indication

NanoBone® bone grafting material is for the reconstruction of bone defects and for augmentation.

Properties

NanoBone® bone grafting material is made up of nanocrystalline hydroxyapatite (HA) which is embedded in a silica gel matrix (amorphous SiO₂). To control the speed of biodegradation, various products with different HA/SiO₂ ratios are offered (retardation of biodegradation with reduction in proportion of SiO₂). The transmission of infections or diseases of animal origin is ruled out as a result of it being manufactured entirely synthetically. NanoBone® products are available as granulate or block, to meet the demands of the application. The porosity of the bone grafting material is approx. 80 %.

Application

NanoBone® | granulate and NanoBone® | block may only be used by qualified professionals with experience in the required surgical techniques. The principles of sterile work and patient medication valid for surgical procedures must be observed.

Preparation:

- Fragments of bone and soft tissue in the defect must be removed from the bone defect to be treated.
- Depending on the nature and location of the bone defect, additional osteosynthetic measures to stabilise the bone grafting material are required.
- NanoBone® bone grafting material may not be implanted into the bone defect dry. For this it is recommended that the granulate or block is mixed or saturated with the patient's own venous blood before implantation. Mixing or saturating with physiological saline solution is also possible.
- When contouring the NanoBone® | block and/or bone bedding, the greatest possible contact surface between the bone grafting material and bone must be produced.

Implantation:

- The area where NanoBone® is implanted must be properly secured with rigid fixation.
- Maximize the contact between the patient's bone and NanoBone® to ensure bone regeneration.
- The bone defect is to be completely filled with NanoBone® | granulate. Avoid overfilling as well as compaction of the bone grafting material in the bone defect.
- When using the NanoBone® | block, contact between the bone bedding and bone grafting material must be guaranteed using additional measures. One option is to secure the block using bone screws and bone plates. The block is pre-drilled for this purpose.
- The filled bone defect must be closed with a primary wound closure ensuring that the defect is covered by the periosteum.

Information

Immobilization of the bone grafting material is necessary for successful bone regeneration. Micro movements can impair bone regeneration and must be avoided using suitable surgical techniques.

As a result of its high porosity, NanoBone® bone grafting material can look like spongy bone in x-rays and therefore immediately after implantation show a corresponding x-ray contrast.

NanoBone® bone grafting material may no longer be used after its use-by date has expired.

Warnings

- There is no data available on the use of NanoBone® bone grafting material with pregnant women and children.
- To avoid contamination, the sterile bone grafting material may only be removed from the packaging immediately before use and is to be applied in a sterile working environment.
- Preparation and resterilisation of unused NanoBone® bone grafting material is not possible, as the destruction of the nanostructure or changes to properties cannot be ruled out during preparation.
- Single use only. Do not attempt to reuse and/or re-sterilize unused NanoBone®.

Interactions

Interactions between NanoBone® bone grafting material and other medicines and drugs are not known.

Side-effects

Side-effects as a result of using NanoBone® bone grafting material are not known.

Contraindications

- Acute and chronic infections in the area of the surgery (soft tissue infections; inflammatory, bacterial bone disease; osteomyelitis); with antibiotic treatment the user is to decide on the use of NanoBone® | granulate by means of a risk-benefits assessment
- Uncontrolled metabolic diseases (e.g. diabetes mellitus)
- Systemic diseases that negatively impact bone regeneration
- Use of medication known to negatively impact bone regeneration
- Immunosuppressive treatment
- Malignant tumours
- Application in areas of open epiphysal plates
- Application not in mechanically instable defects (after stabilization of the defect, e.g. with an osteosynthesis plate the application of NanoBone® is possible),
- In the radiated area and/or before a planned radiation
- Bone defects with no or poor vascularization

Sterility

The bone grafting material is sterilised by irradiation. The expiry date and integrity of the sterile packaging must be checked before preparation.

A copy of the summary of safety and clinical performance (SSCP) can be requested at info@germany@biocomposites.com.

NanoBone® | granulate

CE 0482

NanoBone® | block
Synthetic, biodegradable bone grafting material

Indication

NanoBone® bone grafting material is for the reconstruction of bone defects and for augmentation.

Properties

NanoBone® bone grafting material is made up of nanocrystalline hydroxyapatite (HA) which is embedded in a silica gel matrix (amorphous SiO₂). To control the speed of biodegradation, various products with different HA/SiO₂ ratios are offered (retardation of biodegradation with reduction in proportion of SiO₂). The transmission of infections or diseases of animal origin is ruled out as a result of it being manufactured entirely synthetically. NanoBone® products are available as granulate or block, to meet the demands of the application. The porosity of the bone grafting material is approx. 80 %.

Application

NanoBone® | granulate and NanoBone® | block may only be used by qualified professionals with experience in the required surgical techniques. The principles of sterile work and patient medication valid for surgical procedures must be observed.

Preparation:

- Fragments of bone and soft tissue in the defect must be removed from the bone defect to be treated.
- Depending on the nature and location of the bone defect, additional osteosynthetic measures to stabilise the bone grafting material are required.
- NanoBone® bone grafting material may not be implanted into the bone defect dry. For this it is recommended that the granulate or block is mixed or saturated with the patient's own venous blood before implantation. Mixing or saturating with physiological saline solution is also possible.
- When contouring the NanoBone® | block and/or bone bedding, the greatest possible contact surface between the bone grafting material and bone must be produced.

Implantation:

- The area where NanoBone® is implanted must be properly secured with rigid fixation.
- Maximize the contact between the patient's bone and NanoBone® to ensure bone regeneration.
- The bone defect is to be completely filled with NanoBone® | granulate. Avoid overfilling as well as compaction of the bone grafting material in the bone defect.
- When using the NanoBone® | block, contact between the bone bedding and bone grafting material must be guaranteed using additional measures. One option is to secure the block using bone screws and bone plates. The block is pre-drilled for this purpose.
- The filled bone defect must be closed with a primary wound closure ensuring that the defect is covered by the periosteum.

Information

Immobilization of the bone grafting material is necessary for successful bone regeneration. Micro movements can impair bone regeneration and must be avoided using suitable surgical techniques.

As a result of its high porosity, NanoBone® bone grafting material can look like spongy bone in x-rays and therefore immediately after implantation show a corresponding x-ray contrast.

NanoBone® bone grafting material may no longer be used after its use-by date has expired.

Warnings

- There is no data available on the use of NanoBone® bone grafting material with pregnant women and children.
- To avoid contamination, the sterile bone grafting material may only be removed from the packaging immediately before use and is to be applied in a sterile working environment.
- Preparation and resterilisation of unused NanoBone® bone grafting material is not possible, as the destruction of the nanostructure or changes to properties cannot be ruled out during preparation.
- Single use only. Do not attempt to reuse and/or re-sterilize unused NanoBone®.

Interactions

Interactions between NanoBone® bone grafting material and other medicines and drugs are not known.

Side-effects

Side-effects as a result of using NanoBone® bone grafting material are not known.

Contraindications

- Acute and chronic infections in the area of the surgery (soft tissue infections; inflammatory, bacterial bone disease; osteomyelitis); with antibiotic treatment the user is to decide on the use of NanoBone® | granulate by means of a risk-benefits assessment
- Uncontrolled metabolic diseases (e.g. diabetes mellitus)
- Systemic diseases that negatively impact bone regeneration
- Use of medication known to negatively impact bone regeneration
- Immunosuppressive treatment
- Malignant tumours
- Application in areas of open epiphysal plates
- Application not in mechanically instable defects (after stabilization of the defect, e.g. with an osteosynthesis plate the application of NanoBone® is possible),
- In the radiated area and/or before a planned radiation
- Bone defects with no or poor vascularization

Sterility

The bone grafting material is sterilised by irradiation. The expiry date and integrity of the sterile packaging must be checked before preparation.

A copy of the summary of safety and clinical performance (SSCP) can be requested at info@germany@biocomposites.com.

NanoBone® | granulate

Synthetic, biodegradable bone void filler

Indication

NanoBone® synthetic bone void filler is intended to be used for orthopaedic applications as a filler for gaps and voids that are not intrinsic to the stability of the bony structure. NanoBone® is also intended to fill, augment, and / or reconstruct maxillofacial osseous bone defects, including periodontal, oral, and craniomaxillofacial applications. NanoBone® is indicated to be placed gently into bony voids or gaps of the skeletal system, i.e., extremities, pelvis, and spine, including use in spinal fusion procedures with appropriate stabilising hardware. These defects may be surgically created osseous defects or osseous defects created from traumatic injury to the bone. The product provides a synthetic bone void filler that resorbs and is replaced by bone during the healing process. This medical device is not intended to be used where high tensile, compressive, torsion and/or bending strength are required except when the surroundings have been properly secured mechanically with rigid fixation hardware.

Properties

NanoBone® bone void filler is made up of nanocrystalline hydroxylapatite (HA) which is embedded in a silica gel matrix (amorphous SiO₂). The transmission of infections or diseases of animal origin is ruled out as a result of it being manufactured entirely synthetically. The porosity of the bone void filler is approx. 80 %.

Application

NanoBone® | granulate may only be used by qualified professionals with experience in the required surgical techniques. The principles of sterile work and patient medication valid for surgical procedures must be observed.

Preparation:

- Fragments of bone and soft tissue in the defect must be removed from the bone defect to be treated.
- Depending on the nature and location of the bone defect, additional osteosynthetic measures to stabilise the bone void filler are required.
- NanoBone® bone void filler may not be implanted into the bone defect dry. For this it is recommended that the granulate is mixed or saturated with the patient's own venous blood before implantation. Mixing or saturating with physiological saline solution is also possible.

Implantation:

- The area where NanoBone® is implanted must be properly secured with rigid fixation.
- Maximise the contact between the patient's bone and NanoBone® to ensure bone regeneration.
- The bone defect is to be completely filled with NanoBone® | granulate. Avoid overfilling as well as compaction of the bone void filler in the bone defect.
- The filled bone defect must be closed with a primary wound closure ensuring that the defect is covered by the periosteum.

Information

Immobilisation of the bone void filler is necessary for successful bone regeneration. Micro movements can impair bone regeneration and must be avoided using suitable surgical techniques. As a result of its high porosity, NanoBone® bone void filler can look like spongy bone in x-rays and therefore immediately after implantation show a corresponding x-ray contrast. NanoBone® bone void filler may no longer be used after its use-by date has expired.

Warnings

- NanoBone® is not intended for load-bearing uses.
- There is no data available on the use of NanoBone® bone void filler with pregnant women and children.
- To avoid contamination, the sterile bone void filler may only be removed from the packaging immediately before use and is to be applied in a sterile working environment.
- Preparation and resterilisation of unused NanoBone® bone void filler is not possible, as the destruction of the nanostructure or changes to properties cannot be ruled out during preparation.
- Single use only. Do not attempt to reuse and/or re-sterilise unused NanoBone®.

Interactions

Interactions between NanoBone® bone void filler and other medicines and drugs are not known.

Side-effects

Side-effects as a result of using NanoBone® bone void filler are not known.

Contraindications

- Acute and chronic infections in the area of the surgery (soft tissue infections; inflammatory, bacterial bone disease; osteomyelitis); with antibiotic treatment the user is to decide on the use of NanoBone® | granulate by means of a risk-benefits assessment
- Serious metabolic diseases such as e.g. serious, uncontrollable or poorly controlled diabetes mellitus
- Calcium metabolism disorders
- Malignant tumour
- Treatment with steroids and other pharmaceuticals which interfere with calcium metabolism
- Immunosuppressive treatment
- Endocrinological bone diseases
- Application in areas of open epiphyseal plates
- In the radiated area and/or before a planned radiation
- Bone defects with no or poor vascularisation

Sterility

The bone void filler is sterilised by irradiation. The expiry date and integrity of the sterile packaging must be checked before preparation.

A copy of the summary of safety and clinical performance (SSCP) can be requested at info@germany@biocomposites.com.

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**Biocomposites GmbH**
Biocomposites Ltd
Biocomposites Inc
My Biologics PTY LTD