NanoBone® QD

Bone Graft Substitute

Product Description

NanoBone® QD contains NanoBone® SBX Putty, a synthetic, biodegradable bone grafting material that is composed of non-sintered nanocrystalline hydroxyapatite (HA) which is embedded in a silica gel matrix (amorphous SiO₂). The ratio of HA / SiO₂ has been chosen to optimize the rate of biodegradation for treatment of osseous defects in human bone. NanoBone® QD does not contain any components of animal or human origin, eliminating the possibility of transmission of infection or disease.

NanoBone® QD is supplied in a sterile polymer cartridge with a separate sterile plunger. The cartridge contains NanoBone® granules (0.6-2.0 mm) mixed with a silica gel matrix.

NanoBone® QD does not set in-situ following implantation.

Indications for Use

NanoBone® QD Putty is an implant intended to fill bony voids or gaps of the skeletal system (i.e., extremities, posterolateral spine and pelvis). These osseous defects are surgically created or the result of traumatic injury to the bone and are not intrinsic to the stability of the bony structure. NanoBone® QD Putty resorbs and is replaced with bone during the healing process.

Contraindications

NanoBone® QD is not designed or sold for any use except as indicated. Do not use NanoBone® QD in the presence of any contraindication. NanoBone® QD is contraindicated where the device is intended as structural support in the skeletal system.

Other conditions for which Nano-Bone® QD is contraindicated include:

- Severe vascular or neurological disease
- Uncontrolled diabetes
- Severe degenerative disease
- Severely impaired renal function
- Hypercalcemia, abnormal calcium metabolism
- Existing acute or chronic infections, especially at the site of the operation
- Inflammatory bone disease such as osteomyelitis
- Malignant tumors
- Uncooperative patients who cannot or will not follow post-operative instruction, including individuals who abuse drugs and/or alcohol

Application

NanoBone® QD may only be used by qualified professionals with experience in the required surgical techniques. The principles of sterile work and patient medication valid for surgical procedures must be observed. NanoBone® QD is not intended for load-bearing applications. The area where NanoBone® QD is implanted must be properly secured with rigid fixation. Maximize the contact between the patient's bone and NanoBone® QD to ensure bone regeneration. NanoBone® QD does not set in-situ following implantation. The bone defect must be completely filled with NanoBone® QD. Avoid overfilling as well as compaction of the bone grafting material in the bone defect. The filled bone defect should be closed with a primary wound closure ensuring that the defect is covered by the periosteum. NanoBone® QD may appear radiopaque in x-rays, similar to cancellous bone.

Do not overfill or attempt to pressurize the bony defect site, as this may lead to extrusion of the product beyond the site of its intended application and damage to the surrounding tissues, or may lead to fat embolization or embolization of the product into the bloodstream.

Instructions for Use

1. Open outer (non-sterile) and inner (sterile) packaging.
2. Remove the two wing nuts by turning counterclockwise (Fig. 1).
3. Insert the plunger into the cartridge (Fig. 2).
4. Push the biomaterial out of the cartridge into the bone void (Fig. 3).
5. NanoBone® QD is designed to be used alone.
6. Secure the surgical site after implantation to prevent micromotion and/or implant migration. If NanoBone® QD material is not positioned properly, remove the material and begin again with a new package of material.

Warnings

NanoBone® QD is not intended for load-bearing uses. There is no data available on the use of NanoBone® QD with pregnant women or children.

Single use only. Do not attempt to reuse and/or resterilize unused NanoBone® QD.

Sterility


Contents are STERILE by exposure to gamma radiation unless opened or damaged. Do not use if sterile packaging has been damaged. Do not use if expiration date has been exceeded. Do not resterilize.

Storage, Shelf Life and Disposal

Product should be stored between 5-32°C (40-90°F). The expiration date is printed on the label. DO NOT USE NanoBone® QD AFTER THE EXPIRATION DATE. NanoBone® QD is environmentally friendly. No special disposal is necessary.

Caution

U.S. Federal Law restricts this device to sale by or on the order of a physician. Responsibility for proper selection of patients, for adequate training, for experience in the selection and placement of NanoBone® QD and for the selection of post-operative protocols rests entirely with the physician.

Glossary

| Symbol | Description | Source |
|--------|--|---|
| | Manufactured by. | AAMI/ANSI/ISO 15223-1: 2016, section 5, reference no. 5.1.1 |
| | Contents are STERILE by exposure to gamma radiation. | AAMI/ANSI/ISO 15223-1: 2016, section 5, reference no. 5.2.4 |
| | Do not use if package is damaged. | AAMI/ANSI/ISO 15223-1: 2016, section 5, reference no. 5.2.8 |
| | Do not resterilize. | AAMI/ANSI/ISO 15223-1: 2016, section 5, reference no. 5.2.6 |
| | Single use. Do not reuse. | AAMI/ANSI/ISO 15223-1: 2016, section 5, reference no. 5.4.2 |
| | Store between 40-90 °F. | AAMI/ANSI/ISO 15223-1: 2016, section 5, reference no. 5.3.7 |
| | Consult instructions for use. | AAMI/ANSI/ISO 15223-1: 2016, section 5, reference no. 5.4.3 |
| | CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician. | 21 U.S. Code § 353, paragraph (b)(4)(A) |

Biocomposites GmbH
Fischerweg 421 | 18069 Rostock | Germany
Tel +49 (0) 381 806 99 40
Fax +49 (0) 381 806 99 49 9
Email: infogermany@biocomposites.com
Web: www.biocomposites.com

NanoBone® QD

Synthetic, biodegradable bone grafting material

Indication

NanoBone® QD is for the reconstruction of bone defects which are not essential for the stability of the bone structure and for augmentation. NanoBone® QD can be used for all fields of application for which autologous cancellous bone could be used.

Properties

NanoBone® QD consists of NanoBone® | granulate S39 (0.6-2.0 mm) and a water soluble carrier material containing poloxamer. In NanoBone® | granulate S39, nanocrystalline hydroxyl apatite (HA) is embedded in a silica gel matrix (amorphous SiO₂). The porosity of the granulate is approx. 80 %. The transmission of infections or diseases of animal or human origin is ruled out as a result of it being manufactured entirely synthetically. NanoBone® QD does not harden after implantation.

Application

NanoBone® QD may only be used by qualified professionals with experience in the required surgical techniques. The principles of sterile work and patient medication valid for surgical procedures must be observed.

Preparation:

Fragments of bone and soft tissue in the defect must be removed from the bone defect to be treated.

Depending on the nature and location of the bone defect, additional osteosynthetic measures to stabilise the bone grafting material may be required. Assemble the applicator as shown in Figure 1 "Preparing the applicator".

Implantation:

NanoBone® QD can be used exclusively or in any mixing ratio with autologous cancellous bone.

The area where NanoBone® QD is implanted must be properly secured with rigid fixation.

Maximize the contact between the patient's bone and NanoBone® QD to ensure bone regeneration.

NanoBone® does not set in-situ following implantation.

The bone defect is to be completely filled with NanoBone® QD. Avoid overfilling as well as compaction of the bone grafting material in the bone defect.

The filled bone defect must be closed with a primary wound closure ensuring that the defect is covered by the periosteum.

Information

For successful bone regeneration and in order to avoid the dislocation of NanoBone® QD, a suitable bone position (cavity) must be provided, and the immobilisation of the bone grafting material is necessary. Micro movements can impair bone regeneration and must be avoided using suitable surgical techniques (e.g. using membranes). As a result of its high porosity, NanoBone® bone grafting material can look like spongy bone in x-rays and therefore immediately after implantation show a corresponding x-ray contrast.

NanoBone® QD may no longer be used after its use-by date has expired.

Warnings

NanoBone® QD is not suitable for load-bearing defects. There is no data available on the use of NanoBone® QD with pregnant women and children.

To avoid contamination, the sterile bone grafting material may only be removed from the packaging immediately before use and is to be applied in a sterile working environment.

Preparation of unused NanoBone® QD is not possible as the destruction of the nanostructure or changes to properties cannot be ruled out during preparation.

Single use only. Do not attempt to reuse and/or re-sterilize unused NanoBone® QD.

Interactions

Interactions between NanoBone® QD and other medicines and drugs are not known.

Side-effects

Side-effects as a result of using NanoBone® QD are not known.

Contraindications

- Acute and chronic infections in the area of the surgery (soft tissue infections, inflammatory, bacterial bone diseases, osteomyelitis), with antibiotic treatment the user is to decide on the use of NanoBone® QD by means of a risk-benefits assessment
- Uncontrolled metabolic diseases (e.g. diabetes mellitus)
- Systemic diseases that negatively impact bone regeneration
- Use of medication known to negatively impact bone regeneration
- Immunosuppressive treatment
- Malignant tumours
- Application in areas of open epiphysial plates
- Application not in mechanically instable defects (after stabilization of the defect, e.g. with an osteosynthesis plate the application of NanoBone® is possible),
- In the radiated area and/or before a planned radiation
- Bone defects with no or poor vascularisation

Type and content of the packaging

NanoBone® QD is provided in a sterile applicator to meet the demands of the application. This is packaged in a sterile aluminium protection bag (peel-off packaging). This inner bag is packed in an aluminium protection bag (peel-off packaging). NanoBone® QD must no longer be stored once the inner aluminium protection bag has been opened.

Sterility

NanoBone® QD is sterilised using gamma radiation. Before using NanoBone® QD, its expiration date and the impeccable condition of the sterile packaging must be checked.

A copy of the summary of safety and clinical performance (SSCP) can be requested at infogermany@biocomposites.com.

 Biocomposites GmbH
Fischerweg 421 | 18069 Rostock | Germany
Tel +49 (0) 381 806 99 40
Fax +49 (0) 381 806 99 49 9
Email: infogermany@biocomposites.com
Web: www.biocomposites.com


Figure 1: Preparation of the applicator

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