

SYNICEM™

COMPLETE CONTROL
FROM MIXING TO DELIVERY



Precision mixing & delivery of SYNICEM bone cement

Designed to support your use of SYNICEM low-viscosity bone cements, The SYNIMIX **V** and SERINGUE SYNICEM **GUN** mixing and application systems enhance the mechanical properties of the cement and enable precise application in restricted or hard-to-reach areas. ^{1,2}

SYNIMIX **V**

- ✓ Vacuum mixing eliminates vapours and reduces air bubbles for safer handling and stronger cement
- ✓ Ability to cut cannula to desired length
- ✓ Single-use system with reusable accessory gun
- ✓ Highest standards of quality assurance

SERINGUE SYNICEM **GUN**

- ✓ Designed for quick and controlled injection of bone cement into the femoral cavity
- ✓ Includes mixing bowl and spatula
- ✓ Short or long cannula option
- ✓ Single-use system with reusable accessory gun
- ✓ Highest standards of quality assurance

APPLICATION SYSTEMS

SYNIMIX V

Vacuum mixer and application system for low-viscosity orthopaedic cement



SYNIMED V4

Reusable injection gun for use exclusively with the SYNIMIX V system



SERINGUE SYNICEM GUN

Bone cement application system to inject low-viscosity orthopedic cements



SYNIMED GUN

Reusable gun for use with Seringue SYNICEM GUN



PRODUCT	CONTENT	ORDER CODE
SYNIMIX V	1 mixing body - 1 body cap with mixing rod - 1 plastic tube (3 mm inner diameter and 2 m long) - 1 vapour-absorbing cartridge - 1 funnel - 1 raccord coupling for tubes (3 mm to 6 mm diameter) - 1 pressuriser	880439
SYNIMED V4	Sterilisation box and injection gun	880582
SERINGUE SYNICEM GUN	Syringe body with plunger - syringe lid/cap threaded to fix the cannula - long cannula - short cannula - mixing container - spatula - funnel	880447
SYNIMED GUN	Sterilisation box and injection gun	880593

Find out more at biocomposites.com

References: 1. Vacuum mixer and application system for bone cement, Instructions for Use. 2. Application system for surgical cement, Instructions for Use.

©2025, Biocomposites is a trademark/registered trademark of Biocomposites Ltd. SYNICEM is a trademark/registered trademark of Biocomposites International Ltd. All rights reserved. No unauthorised copying, reproduction, distributing, or republication is allowed unless prior written permission is granted by the owner, Biocomposites Ltd/Biocomposites International Ltd.

SYNIMIX V, SERINGUE SYNICEM GUN, SYNIMED V4 and SYNIMED GUN are manufactured by Synergie Ingénierie Médicale S.A.R.L.

Synergie Ingénierie Médicale S.A.R.L. (Synimed), Zone Artisanale de l'Angle, 19370 Chamberet, France.