

SYNICEM CLOUS

ANTIBIOTIC COATED INTRAMEDULLARY NAIL

The SYNICEM **CLOUS** is an antibiotic coated nailing system that provides stability to infected femoral and tibia fractures. The addition of gentamicin helps to treat the infection.¹

- ✓ Infected diaphyseal fractures
- ✓ Infected diaphyseal pseudoarthrosis
- ✓ Diaphyseal osteomyelitis

Core nail material:
316 LVM Stainless Steel

Outer coating:
3.38% gentamicin sulfate
impregnated PMMA

Coating thickness:
1mm

SYNICEM CLOUS

ANTIBIOTIC COATED INTRAMEDULLARY NAIL

Available in a wide range of diameters and lengths

TIBIAL NAIL

Model	Order Code	Model	Order Code	Model	Order Code
Ø9 x 250mm	882377	Ø10 x 250mm	882388	Ø11 x 250mm	882399
Ø9 x 280mm	882379	Ø10 x 280mm	882390	Ø11 x 280mm	882401
Ø9 x 310mm	882381	Ø10 x 310mm	882392	Ø11 x 310mm	882403
Ø9 x 340mm	882383	Ø10 x 340mm	882394	Ø11 x 340mm	882405
Ø9 x 370mm	882385	Ø10 x 370mm	882396	Ø11 x 370mm	882407


FEMORAL NAIL

Model	Order Code	Model	Order Code	Model	Order Code
Ø10 x 320mm	881551	Ø12 x 320mm	881570	Ø14 x 360mm	881590
Ø10 x 340mm	881552	Ø12 x 340mm	881571	Ø14 x 380mm	881591
Ø10 x 360mm	881553	Ø12 x 360mm	881572	Ø14 x 400mm	881592
Ø10 x 380mm	881554	Ø12 x 380mm	881573		
Ø10 x 400mm	881555	Ø12 x 400mm	881574		

Note: Additional nail diameters and lengths are available subject to minimum order quantities.

SYNICEM CLOUS iSET

The SYNICEM CLOUS iSET instrument set is designed to support a surgeon's use of the SYNICEM CLOUS antibiotic coated intramedullary nail.²



Order Code
882074

References: 1. Synicem Cloud, Instructions for Use. 2. Synicem Cloud iSET, Instructions for Use.

©2025, Biocomposites Ltd. Synimed and Synicem are trademarks/registered trademarks of Synergie Ingénierie Médicale S.A.R.L. All rights reserved. No unauthorised copying, reproduction, distributing, or republication is allowed unless prior written permission is granted by the owner, Biocomposites Ltd.

The Synicem Cloud and Synicem Cloud iSet are manufactured by Synergie Ingénierie Médicale S.A.R.L. (Synimed), Zone Artisanale de l'Angle, 19370 Chamberet, France.

Patents granted: EP 2 594 297 B1, US 9,770,273 B2

MA0520R1